

IN VITRO TRANSCRIPTION

Protocol for the *in vitro* transcription for RNAi libraries. The protocol is optimized for high-throughput *in vitro* transcription.

IVT reactions are performed in 96-well plates from Thermo-Scientific (cat# AB-0796) using a waterbath.

The IVT reaction are performed according to the manufacturer's description depending on the Supplier. Reactions should be performed in a volume of 50 μ l.

We use 15 μ l of the T7 containing 2ndPCR product as template. Incubate for several hours to overnight.

- DNase treatment

	Per 50μl reaction
RNase free water (Acros Organics)	19.4 μ l
DNaseI (Fermentas)	0.6 μ l

- Incubate at room temperature for 30min
- Ethanol precipitation

	Per 50μl reaction
5M NH ₄ OAc (5M, Ambion)	8 μ l
Ethanol absolut (Riedel-DeHaen)	200 μ l

- Incubate at 2-24h at -20°C
- Spin plates for 30min at 4000rpm
- Decant supernatant
- Dry pellet 1,5-2h
- Resuspend pellet in 100 μ l 10mM Tris pH7.0 (Ambion)
- Quality check on a gel (dilute IVT 1:50 and load 5 μ l dsRNA + 15 μ l loading dye)

Distribution: internal OK